

CE

WYLEDKL4

LED Operating Light



User-friendly and Best Value for Money



Osram LED lighting deliver cool, white bright light;



Control panel, switch on/off & 10 grade intensity adjustment




Optional stainless base, 7.8ah lithium battery, last 4 hours once fully charged

WYLEDKL4

Floor Standing LED Operating Light

Premium Quality, Surprising Price




Smooth light focusing system

General Features

- High-intensity of 100,000 lux at 1 meter;
- 4 500K color temperature offers natural and true coloring;
- Accurate color rendition;
- Super deep-cavity illumination;
- Excellent shadow control;
- Energy efficiency;
- Germany imported Osram LED chips – each rated with a 50,000 hour bulb life;
- Extremely smooth focus mechanism with the focusing knob;
- Less maintenance, provide best outcome for hospitals.

Technical Specifications

- Illumination: 40,000-100,000 lux, adjustable
- Color temperature: 4500±500K
- Color rendering index: 95%
- Light field size: 150-300mm
- Illumination depth: 800mm
- Bulb power: 1W
- Diameter of lamp head: 500mm
- Qty of LEDs: 20
- Temperature rise at surgeon's head: Less than 2°C



All the measurement is millimetre.