



▲ XSZ-126III

- ① Sharp Image Display with Excellent Optical Function.
- ② Small Space Occupancy with Integral Stand Design, Low Environment Requirement with Anti-mould Technology.
- ③ Comfortable Operation with Ergonomic Design Principle.

XSZ-126 Series Biological Microscope Outfits

Items	Specification	XSZ-126III
Optical System	Finite Optical System	●
Eyepiece	WF10X/18mm	●
	WF16X	○
Viewing Head	Seidentopf Binocular Viewing Head Inclined at 30° , Interpupillary 48-75mm	○
	Seidentopf Trinocular Viewing Head Inclined at 30° , Interpupillary 48-75mm	●
Objective	Achromatic 4X、10X、40X(S)、100X(S、 OIL)	●
	Semi-Plan 4X、10X、40X(S)、100X(S、 OIL)	○
	Plan 4X、10X、40X(S)、100X(S、 OIL)	○
Nosepiece	Backward Quadruple Nosepiece	●
Stage	Double Layers Mechanical Stage 140X140mm, Moving Range 75X50mm	●
Condenser	NA1.25 Abbe Condenser	●
	Dark Field Condenser (Dry, Oil)	○
Focusing System	Coaxial Coarse and Fine Adjustment, Fine Division 0.002mm, Moving Range 25mm	●
Illumination	Halogen Lamp 6V/20W	○
	S-LED Illumination	●
	S-LED Illumination with Recharger	○
	Plane-concave Mirror	○
Fluorescent Attachment	Fluorescent Illumination, LED lamp	○
	Fluorescent Illumination, 100W HBO lamp	○

Note: " ● "In Table Is Standard outfits, " ○ " Is Optional Accessories.

Packing Size: 390mm×260mm×490mm Gross Weight: 10.0kgs Net Weight: 8.0kgs

宁波盛恒光电有限公司

NINGBO SHENG HENG OPTICS & ELECTRONICS CO.,LTD

⑨ 宁波市海曙区高桥工业区秀丰路 587 号
No.587 Xiu Feng Road Gao Qiao Industrial Estate Haishu District
Ningbo China 315174

☎ +86-574-88046016/88046015/88054001

☎ +86-574-88447692/88446043

🌐 www.microscopechina.com

✉ navite@navite.com sales@navite.com yhj@navite.com

